

# Cytoflex 2-L Cheat Sheet

Excitation	Filters	Fluorochromes
488nm	525/40 BP	FITC, Alexa488, CFSE, GFP
	585/42 BP	PE, PI
	610/20 BP	ECD, PE-Texas Red, PE-CF594, PI
	690/50 BP	PC5, PC5.5, PerCP, PerCP-CY5.5
	780/60 BP	PE-CY7, PC7
637nm	660/20 BP	APC, Alexa Fluor 647
	712/25 BP	APC-Alexa700, Alexa Fluor 700
	780/60 BP	APC-Alexa750, APC-H7, APC-CY7, APC-Efluor780

- Blue Laser: 50mW
- Red Laser: 50mW

Computer based training:

<https://www.beckman.com/service/training/on-demand/flow-cytometry/cytoflex-online-curriculum>

Fluorescence Spectrum Analyzer:

<https://www.beckman.com/flow-cytometry/fluorescence-spectrum-analyzer>

Maintenance how to video:

<https://www.mybeckman.ca/resources/videos/products/cytoflex-special-procedure>

For training or assistance

Contact : Mike 412 475 6312  
[meyere@upmc.edu](mailto:meyere@upmc.edu)